

NEW!!

Diamond Soft Twist

Teklad 7177C

A combination of Diamond Soft bedding with Diamond Twist enrichment!



PRODUCT DESCRIPTION

Diamond Soft Twist are the newest bedding and enrichment product from Teklad. This new product combines our popular Diamond Soft bedding plus ½ inch pieces of Diamond Twists.

This contact bedding is the perfect combination of specialty processed hardwood paper pulp making this a highly absorbent paper product. Due to its white coloration, it enables the observation of any discharges from the animal after birth or surgery.

To enhance animal welfare benefit and provide additional enrichment, we incorporated our popular Diamond Twist enrichment product. These ½ inch pieces of thick rolled tissue paper offer low dust, ensuring a clean and healthy environment. This product is shreddable and does not wind around rodents legs or pups. The paper material will not clog drains or disposals.

PRODUCT FEATURES AND BENEFITS

- Utilizes the softest contact bedding material to ensure animal comfort
- Nestable design promoting burrowing and nesting behaviors
- Eliminates the need for additional enrichment
- Encourages natural foraging behaviors in rodents
- Bright white color provides improved animal observations
- Highly absorbent material effectively controls ammonia levels
- Available irradiated and autoclavable packaging

PACKAGING

Diamond Soft Twist is available in a 25 lb perforated autoclavable bag.

Find out more

Contact your Inotiv representative
or visit us at [inotiv.com](https://www.inotiv.com)