

TD.00278

Purified Diet, Low Residue

Formula

g/Kg

Whey Protein Isolate	205.0
DL-Methionine	2.2
Dextrose, monohydrate	528.39
Maltodextrin	150.0
Soybean Oil	50.0
Mineral Mix, AIN-93G-MX (94046)	35.0
Potassium Phosphate, monobasic	8.4
Calcium Carbonate	3.0
Vitamin Mix, AIN-93-VX (94047)	15.0
Choline Bitartrate	3.0
TBHQ, antioxidant	0.01

Footnote

Diet uses highly digestible ingredients and excludes a source of fiber. These features are often desired to reduce fecal output following gastrointestinal surgeries.

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	19.3	20.2
Carbohydrate	64.8	67.9
Fat	5.1	11.9
Kcal/g	3.8	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotiv.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Low Residue
- Fiber Free

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotiv.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotiv.com