

# TD.02028

## Atherogenic Rodent Diet

Formula	g/Kg
Casein	195.0
DL-Methionine	3.0
Sucrose	325.46
Corn Starch	150.0
Anhydrous Milkfat	210.0
Cholesterol	12.5
Cholic Acid	5.0
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.04

### Footnote

This is a modification of diet TD.88137 to increase added cholesterol to 1.25% (total in the diet is 1.3%) and add 0.5% cholic acid. TD.96121 is a version without cholic acid. For questions about this diet or modifications [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

### Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	17.3	15.4
Carbohydrate	46.9	41.9
Fat	21.2	42.6

**Kcal/g 4.5**

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

### Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



### Key Features

- Purified Diet
- Modification of TD.88137
- Cholesterol
- Cholic Acid

### Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



### Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

### Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

### Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)