

TD.10885

45% Fat Kcal Diet (0.2% TC)

Formula	g/Kg
Casein	195.0
L-Cystine	3.0
Sucrose	340.0
Corn Starch	55.46
Maltodextrin	60.0
Anhydrous Milkfat	210.0
Cholesterol	1.5
Soybean Oil	20.0
Cellulose	50.0
Mineral Mix, AIN-93G-MX (94046)	43.0
Vitamin Mix, AIN-93-VX (94047)	19.0
Choline Bitartrate	3.0
TBHQ, antioxidant	0.04

Footnote

Designed with similarities to TD.88137 Western Diet with 21% anhydrous milk fat and 34% sucrose. 2% soybean oil is included to supplement essential fatty acids. Modified corn starch (maltodextrin) is added for improved pellet quality. Total cholesterol of 0.2%. For questions about this diet or modifications askanutritionist@inotivco.com

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	17.3	14.8
Carbohydrate	47.5	40.6
Fat	23.2	44.6
Kcal/g	4.7	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only
Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Modification of TD.88137
- Improved Pellet Quality
- Improved Essential Fatty Acids

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com