

Research Models and Services Teklad Bedding and Enrichment

Pelleted Paper Bedding

Teklad 7084

Teklad 7084 Pelleted Paper Bedding is an

environmentally friendly choice in laboratory animal bedding material.

PRODUCT DESCRIPTION

Pelleted Paper Bedding is made from pre-consumer paper, mixed with wood dust as a pelleting aid, so it is very clean, highly absorbent, and non-toxic. Originally designed as a non-contact bedding and cat litter, pelleted paper has been increasingly used as a contact bedding with excellent results, particularly in superior absorbency and odor control. It also will not plug automated watering valves.

Unlike most other pelleted paper bedding products, which are made up mostly of post-consumer recycled newspaper, Teklad 7084 Pelleted Paper Bedding is made from pre-consumer recycled paper, most of which was never printed upon. This is a very important distinction, because this product will not contain any bits of plastic, string, staples, high ink content, or other contaminants that can be introduced in the post-consumer recycling process. Pre-consumer paper is also never subjected to any chemical de-inking solvents. For the small percentage of raw material that we use that does have printing on it, only soy ink is used. The net result is a much safer, cleaner, and lighter colored product that will not turn your white rodents gray.



PRODUCT FEATURES AND BENEFITS

- Environmentally friendly
- Biodegradable and incinerable
- Highly absorbent
- Very low dust
- Non-toxic (No post-consumer contaminants, chemicals or solvents)
- Excellent odor control
- Lighter color

PRODUCT FORMATS

7084 Pelleted Paper Bedding is produced in a plant that is dedicated solely to the production of laboratory animal bedding products. This plant employs stringent quality assurance, housekeeping and sanitation procedures consistent with what the research community demands and expects. Laboratory analysis of the finished product is performed at least annually, and provides information on potential contaminants such as pesticides, heavy metals, yeast and molds. These certificates of analysis are available upon request.

PACKAGING

Teklad 7084 Pelleted Paper Bedding is packaged in 40 lb. (18.17 kg) autoclavable bags.

Find out more

Contact your Inotiv representative or visit us at **inotivco.com**