

# TD.83171

## Vitamin Mix without A, D, E, Choline

Formula	g/Kg
Biotin	0.08
Calcium Pantothenate	13.22
Folic Acid	0.4
Inositol	22.02
Vitamin K <sub>3</sub> , menadione	9.92
Niacin	19.82
p-Aminobenzoic Acid	22.02
Pyridoxine HCl	4.4
Riboflavin	4.4
Thiamin (81%)	4.4
Vitamin B <sub>12</sub> (0.1% in mannitol)	5.94
Vitamin C, ascorbic acid, coated (97.5%)	203.32
Corn Starch	690.06

### Footnote

This modification of general purpose, well-fortified vitamin mix 40060 is designed for use at 0.5% (5.0 g/kg), to permit adjustment of vitamin A, D, E, and choline in the diet. This mix includes a source of vitamin C.

### Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)



### Key Features

- Vitamin Mix
- Vitamins A, D & E
- Choline

### Key Planning Information

- Products are made fresh to order
- Lead time 2 weeks
- Store product at 4°C or lower
- Use within 6 months (applicable to most mixes)
- Box labeled with product name, manufacturing date, and lot number



### Product Specific Information

- Powder
- Minimum order 500 grams

### Options (fees will apply)

- Rush order (pending availability)
- Vacuum packaging (1 and 2 Kg)

### Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

### Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)