TD.88239

Potassium Deficient Diet

Formula	g/Kg
Casein	200.0
DL-Methionine	3.0
Sucrose	493.7586
Corn Starch	150.0
Corn Oil	50.0
Cellulose	50.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.01
Calcium Phosphate, dibasic	20.3
Calcium Carbonate	10.0
Sodium Chloride	7.4
Magnesium Oxide	2.8
Magnesium Sulfate	2.0
Ferric Citrate	0.5
Manganous Carbonate	0.13
Zinc Carbonate	0.06
Cupric Carbonate	0.02
Chromium Potassium Sulfate, dodecahydrate	0.02
Sodium Iodate	0.001
Sodium Selenite, pentahydrate	0.0004

Footnote

Excludes added potassium. Estimated background potassium is 15-30 ppm. Calculated values for some other minerals: 1% calcium, 0.6% phosphorus, 0.2%, magnesium, 0.3% sodium and 0.45% chloride. TD.88238 is a possible ingredient matched control with 1% potassium. For questions about this or other potassium adjusted diets askanutritionist@inotivco.com.

Selected Nutrient Information ¹

		% by weight	% kcal from
Protein		17.7	19.1
Carbohydrate		63.3	68.3
Fat		5.2	12.6
Kcal/g	3.7		

Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

askanutritionist@inotivco.com



Key Features

- Purified Diet
- Potassium
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
 More frequent replacement may be advised
- Lead time:
 - · 2 weeks non-irradiated
 - · 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

teklad@inotivco.com