

TD.90018

40% Protein Diet

Formula

	g/Kg
Casein	460.0
Sucrose	231.82
Corn Starch	200.0
Corn Oil	50.0
Cellulose	15.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.01
Mineral Mix, Ca-P Deficient (79055)	13.37
Calcium Phosphate, dibasic	9.6
Calcium Carbonate	10.2

Footnote

The 40% protein diet in a popular series of adjusted protein diets. These diets are isocaloric and are matched for calcium (0.7%) and phosphorus (0.54%).

Selected Nutrient Information ¹

	% by weight	% kcal from
Protein	40.0	42.6
Carbohydrate	41.6	44.3
Fat	5.5	13.1
Kcal/g	3.8	

¹ Values are calculated from ingredient analysis or manufacturer data

Collaborate with a Nutritionist

- askanutritionist@inotivco.com

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | inotivco.com



Key Features

- Purified Diet
- Protein Adjusted Series
- Isocaloric
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
More frequent replacement may be advised
- Lead time:
 - 2 weeks non-irradiated
 - 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/tekklad-orders
- tekkladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- tekklad@inotivco.com