

TD.92052

Swine Diet (15% Lard, 1.5% Chol.)

Formula	g/Kg
Corn	283.584
Wheat Middlings	250.0
Soybean Meal	222.0
Alfalfa Meal (17%), dehydrated	30.0
Lard	150.0
Cholesterol	15.0
Dicalcium Phosphate, FG (18.5% P, 21% Ca)	14.0
Calcium Carbonate, FG (38%)	11.2
Mineral Mix, Swine (84246)	5.0
DL-Methionine, FG (99%)	1.6
L-Lysine HCl, FG (78%)	0.6
Vitamin Mix, AIN-76A (40077)	15.0
Choline Chloride, FG (60%)	2.0
Vitamin B ₁₂ (0.1% in mannitol)	0.01
Vitamin K, MSB complex	0.006

Footnote

A swine diet consisting of natural ingredients, with 15% lard and 1.5% cholesterol. Similar to TD.93296 (2% cholesterol) and TD.92033 (1.2% cholesterol). Powder is the recommended form; pellets are very soft and crumbly.

Collaborate with a Nutritionist

- askanutritionist@inotivco.com



Key Features

- Natural Ingredient Diet
- Atherogenic
- Lard
- Cholesterol

Key Planning Information

- Products are made fresh to order
- Lead time 2 weeks
- Store product at 4°C or lower
- Use within 6 months (applicable to most mixes)
- Box labeled with product name, manufacturing date, and lot number



Product Specific Information

- Powder
- Minimum order 500 grams

Options (fees will apply)

- Rush order (pending availability)
- Vacuum packaging (1 and 2 Kg)

Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

Price Quotes | Order Status

- teklad@inotivco.com