

# TD.93296

## Swine Diet (15% Lard, 2% Cholesterol)

Formula	g/Kg
Corn	276.684
Wheat Middlings	250.0
Soybean Meal	223.9
Alfalfa Meal (17%), dehydrated	30.0
Lard	150.0
Cholesterol	20.0
Dicalcium Phosphate, FG (18.5% P, 21% Ca)	14.0
Calcium Carbonate, FG (38%)	11.2
Mineral Mix, Swine (84246)	5.0
DL-Methionine, FG (99%)	1.6
L-Lysine HCl, FG (78%)	0.6
Vitamin Mix, AIN-76A (40077)	15.0
Choline Chloride, FG (60%)	2.0
Vitamin B <sub>12</sub> (0.1% in mannitol)	0.01
Vitamin K, MSB complex	0.006

### Footnote

A swine diet consisting of natural ingredients, with 15% lard and 2% cholesterol. Similar to TD.92052 (1.5% cholesterol) and TD.92033 (1.2% cholesterol). Powder is the recommended form; pellets are very soft and crumbly.

### Selected Nutrient Information<sup>1</sup>

	% by weight	% kcal from
Protein	17.2	19.2
Carbohydrate <sup>2</sup>	33.4	37.3
Fat	17.3	43.5
<b>Kcal/g</b>	<b>3.6</b>	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

<sup>2</sup> Estimated digestible carbohydrate

### Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)



### Key Features

- Natural Ingredient Diet
- Atherogenic
- Lard
- Cholesterol

### Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



### Product Specific Information

- Powder
- Minimum order 3 Kg
- Irradiation available upon request

### Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

### Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)

Teklad diets are designed & manufactured for research purposes only

Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)

032823