# TD.94059

## Purina #5015, Cocoa Butter, Cholesterol, etc.

Formula	g/Kg
5015, PMI Mouse Diet	750.0
Casein	75.0
Dextrose, monohydrate	30.0
Sucrose	16.25
Dextrin	16.25
Cocoa Butter	75.0
Cholesterol	12.5
Cellulose	12.5
Mineral Mix, AIN-76 (170915)	8.75
Vitamin Mix, Teklad (40060)	2.5
Choline Chloride	1.25



## **Key Features**

- Mix of Natural & Refined Ingredients
- Atherogenic
- Cholesterol
- Excludes Cholate

## **Key Planning Information**

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
  More frequent replacement may be advised
- Lead time:
  - · 2 weeks non-irradiated
  - · 4 weeks irradiated



# **Product Specific Information**

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

## Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

askanutritionist@inotivco.com

**Footnote** 

# Selected Nutrient Information <sup>1</sup>

		% by weight	% kcal from	
Protein		20.8	21.7	
Carbohydrate		39.1	41.0	
Fat		15.8	37.3	
Kcal/g	3.8			

A modification of TD.88051/TD.90221 (Paigen hybrid atherogenic diet) without sodium

cholate. The diet contains approx. 15.8% fat (half from cocoa butter) and 1.25%

cholesterol. For guestions about this diet or other modifications

<sup>1</sup>Values are calculated from ingredient analysis or manufacturer data

## **Collaborate with a Nutritionist**

askanutritionist@inotivco.com

# Place Your Order (USA & Canada)

Please Choose One

- www.inotivco.com/teklad-orders
- tekladorders@inotivco.com
- 800.483.5523

## **Price Quotes | Order Status**

• teklad@inotivco.com