

# TD.95027

## Low Calcium Diet (0.4% P)

Formula	g/Kg
Casein	200.0
L-Cystine	3.0
Sucrose	342.188
Corn Starch	320.0
Soybean Oil	60.0
Cellulose	40.0
Mineral Mix, Ca-P Deficient (79055)	13.37
Potassium Phosphate, monobasic	11.43
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.012

### Footnote

Casein contributes approximately 0.01% calcium or less. Diet contains about 0.4% phosphorus and 2200 IU vitamin D/kg diet. TD.97191 is a possible ingredient matched control with ~0.6% calcium. For questions about this or other calcium adjusted diets askanutritionist@inotivco.com.

### Selected Nutrient Information <sup>1</sup>

	% by weight	% kcal from
Protein	17.7	18.6
Carbohydrate	63.4	66.7
Fat	6.2	14.7
Kcal/g	3.8	

<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data

### Collaborate with a Nutritionist

- [askanutritionist@inotivco.com](mailto:askanutritionist@inotivco.com)

Teklad diets are designed & manufactured for research purposes only  
Teklad™ Custom Diets, Madison, WI | [inotivco.com](http://inotivco.com)



### Key Features

- Purified Diet
- Calcium
- Phosphorus

### Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week  
*More frequent replacement may be advised*
- Lead time:
  - 2 weeks non-irradiated
  - 4 weeks irradiated



### Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

### Options (fees will apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### Place Your Order (USA & Canada)

Please Choose One

- [www.inotivco.com/teklad-orders](http://www.inotivco.com/teklad-orders)
- [tekladorders@inotivco.com](mailto:tekladorders@inotivco.com)
- 800.483.5523

### Price Quotes | Order Status

- [teklad@inotivco.com](mailto:teklad@inotivco.com)