

# Certified Diamond Dry Bedding

## Teklad 7070C

Teklad 7070C Certified Diamond Dry is a premium, bright white contact bedding material made from pure cellulose paper.

### PRODUCT DESCRIPTION

Highly defined and highly absorbent, the product is bright white in color to aid in observation of discharges. The overall dust content is very low for a refined white paper material. This product's absorbency helps to minimize ammonia levels. The fluffier nature of this paper material allows for significantly more nesting than with traditional wood chip.

### PRODUCT FEATURES AND BENEFITS

- Highly-defined raw material
- Premium white color
- Ease of observation for health monitoring
- Superior comfort and absorbency
- Promotes nesting
- Ammonia control
- Dispenses easily, including via bulk and automated handling systems
- Monitored for quality and consistency throughout production
- Produced in a plant dedicated to production of lab animal bedding
- Certified and available in both autoclavable and irradiated format



### PRODUCT FORMATS

Teklad 7070C Certified Diamond Dry is monitored for quality and consistency throughout production, and the material is certified for pesticides and heavy metals. Certificates of analysis available on a per-lot basis. It is also available irradiated (7970C).

### PACKAGING

- Teklad 7070C Certified Diamond Dry Bedding packaged in 20 lb. (9.07 kg) autoclavable bags.
- Teklad 7970C Certified Irradiated Diamond Dry Bedding packaged in 20 lb. (9.07kg) cap sacs; the cap sac has a heat-sealed plastic liner within a multi-wall paper bag to ensure purity and integrity post irradiation.

**Find out more**  
Contact your Inotiv representative  
or visit us at [inotivco.com](https://www.inotivco.com)